

# DispozABLE™ Spacer

Low cost, fast response spacer  
designed for emergency use.

- + Disposable design
- + For emergency use - ambulances, schools or GP practices
- + On device pictorial instructions
- + Easy to use



## Essential in emergency situations

Designed for ambulances, schools or GP Practices.

## Disposable and convenient design for single use

Can be used for up to 7 days if required to follow a patient through their treatment journey.

## Simple to assemble

The lid and base are the same, and therefore interchangeable.

## Easy to use

There are pictorial instructions on each individual device, making it simple to use.

## Spirometry reversibility

The DispozABLE Spacer is widely used for Spirometry reversibility.

# DispozABLE™ Spacer

DispozABLE™ Spacer is a low cost, fast response spacer for delivering medication. Designed to be used with ease by non-professionals in an emergency situation.

DispozABLE™ Spacer is a single use spacer; its simple intuitive design with interchangeable lid / base allows for quick assembly. Clear pictorial instructions are displayed on device for immediate identification, adding to the element of ease.

## Cleaning and Maintenance

DispozABLE Spacer is a single use device and can be disposed of in general household waste.

## Technical Specifications

	DEVICE
Dimensions (mm)	ø84 x 107
Weight (g)	15
Volume (ml)	227
Resistance to flow	Free flow

## Performance Data

Fine particle dose (salbutamol)	Delay 0 seconds = 54.5cg (SD 4.8)
	Delay 1 second = 37.7cg (SD 7.6)

PART	DESCRIPTION	QTY
3609550	DispozABLE Spacer	10
3609551	DispozABLE Spacer	200

## Clement Clarke International Ltd

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK  
+44 (0)1443 804999 • [resp@clement-clarke.com](mailto:resp@clement-clarke.com) • [www.clement-clarke.com](http://www.clement-clarke.com)

